



ADPS™

JAK2 V617F Mutation Test kit (RUO)

Detection of G1849T mutation in codon 617 of the JAK2 Gene
Instruction for Use



Background

Myeloproliferative neoplasms (MPN) are group of rare blood cancers derived from myeloid stem cells and increased risk for thrombosis and secondary leukemic transformation. Three major diseases constitute MPN: polycythemia vera (PV), essential thrombocythemia (ET) and primary myelofibrosis (PMF), which are also characterized by their unique somatic mutations

Recognizing the special diagnostic challenge that these diseases represent, a group of pathologists and clinicians sponsored by the World Health Organization (WHO) created the MPN category to provide a less restrictive view of myeloid disorders, which in some instances clearly overlap. The WHO group proposed that the new MPN category would allow for more focused clinical and laboratory investigations of myeloid proliferation, abnormal proliferation, and dysplasia.

Janus kinase 2 (JAK2) is a member of a group of cytoplasmic tyrosine kinases are that involved in the transduction of signals from growth factor receptors.¹ JAK2 is a key player in the JAK–signal transducer and activator of transcription (STAT) pathway. On autophosphorylation following activation via ligand binding, JAK2 recruits STAT molecules, which are then phosphorylated and translocate to the nucleus to act as transcription factors. This mutation disrupts the autoinhibitory JH2 domain of JAK2, which leads to the constitutive activation of JAK2 and, subsequently, uncontrolled proliferation.

A somatic gain-of-function mutation of the JAK2 kinase gene on chromosome 9 has been identified in more than 50% of individuals with MPN. The mutation is an invariant G-to-T transversion in exon 14 resulting in a valine-to-phenylalanine substitution in codon 617 (JAK2V617F). JAK2 V617F mutation (G1849T) is present in a number of myeloproliferative Neoplasms such as PV (~95%), ET (~60%), PMF (~50%).

The incidence of the V617F mutation in patients with bcr-abl– myeloproliferative disorders varies, but defining the presence or absence of this mutation is now part of clinical diagnostic algorithms.

JAK inhibitors that can treat diseases caused by mutations of JAK2 V617F have been evaluated for efficacy and safety through several clinical trials. It has entered the phase trial and is recognized for its excellent efficacy and stability. The Jakafi (Ruxolitinib Phosphate) is approved for use in adults to treat Myelofibrosis (Primary myelofibrosis, Post-polycythemia vera myelofibrosis, Post-essential thrombocythemia myelofibrosis) and Polycythemia vera in patients who cannot be treated with or have not gotten better with hydroxyurea.

However, our understanding of an individual patient's cancer is often limited by tumor accessibility because of the high risk and invasive nature of current tissue biopsy procedures. “Liquid biopsy”, the analysis of circulating tumor DNA (ctDNA), offers a new source of cancer-derived materials that may reflect the status of the disease better and thereby contribute to more personalized treatment. In addition, it is critical to develop more sensitive technologies for detection of mutation.

The ADPS™ JAK2 V617F Mutation Test kit (RUO) is a qualitative PCR assay designed to detect V617F substitution mutation in codon 617. This kit can detect up to 0.01% of the JAK2 V617F mutation in peripheral blood samples or whole blood sample and enables qualitative detection.

Intended Use

The ADPS™ JAK2 V617F Mutation Test kit (RUO) is a real-time PCR assay for qualitative detection of somatic mutation (G1849T) in codon 617 of the JAK2 gene. Genomic DNA extracted from whole blood or circulating free DNA extracted from liquid biopsy could be used. The kit is for Research Use Only (RUO), and intended to be used by trained professionals in a laboratory environment.

Principles of the Procedure

The kit adopts Allele-specific PCR (AS-PCR) and real-time PCR technologies, which comprises specific primers and fluorescent probes to detect *JAK2 V617F* mutation in human plasma DNA samples. During the nucleic acid amplification, the targeted mutant DNA is matched with the bases at 3' end of the primer, amplified selectively and efficiently, then the mutant amplicon is detected by fluorescent probes labeled with FAM. While the wild-type DNA cannot be matched with specific primers, there is no amplification occurs.

The kit is composed of **JAK2 V617F Master Mixture**, **ADPS™ smart DNA polymerase**, **JAK2 V617F Positive Control** and **nuclease-free water**.

- 1) The contents in **JAK2 V617F Master Mixture** formed a mutation detection system and an Internal control system. The mutation detection system includes primers and FAM-labeled probe specific for designated **JAK2 V617F** mutation. The internal control system contains the primers and Quasar 670-labeled probe for a **JAK2** region of genomic DNA without known mutations and polymorphism, to check the presence of inhibitors and monitor the accuracy of experimental operation.
- 2) The **JAK2 V617F Positive Control** contains recombinant plasmid DNA that carries fragments of the **JAK2** gene with V617F mutation.
- 3) The **ADPS™ smart DNA polymerase** contains ADPS DNA polymerase for PCR amplification

Kit Contents

This kit contains the following materials (see Table 1):

Table 1. Kit Contents (for 24 tests)

Contents	Main Ingredient	Quantity
JAK2 V617F Master Mixture	Primers, Probes, dNTPs, Primers/probe for IC*	240 µL/tube ×1
ADPS™ smart DNA polymerase	ADPS™ DNA Polymerase (1 U/µL)	12 µL/tube ×1
JAK2 V617F Positive Control	Recombinant plasmid (JAK2 V617F) (3x10 ⁴ copies/5 µL)	60 µL/tube ×1
Nuclease-free Water	PCR grade water (for no template control, NTC)	1500 µL/tube ×1

*IC: Internal control for PCR

The **JAK2 V617F Master Mixture** contains all reagents except **ADPS™ smart DNA polymerase** that is supplied in a separate tube and has to be added before the assay. The detailed detection information is listed in Table 2.

Table 2. Detection Information

Reagent	Mutation detected	Fluorescent Signal	
		FAM	Quasar 670**
JAK2 V617F Master Mixture	V617F	V617F	IC*

*IC: Internal control for PCR

**Alternative dye: Quasar 670 - CY5

Storage and Stability

The kit requires shipment on frozen ice packs. All contents of the kit should be stored immediately upon receipt at -20±5 °C and protected from light.

The shelf-life of the kit is 18 months, and 90 days of opened reagents for the kit. **JAK2 V617F Master Mixture** and reaction Master Mixtures should be protected from prolonged exposure to light. Reaction Master Mixtures must be stored at 2°C to 8°C in the dark. The prepared samples and controls must be added within 1 hour of preparation of the reaction Master Mixtures. Amplification must be started within 1 hour from the time that the processed samples and controls are added to the reaction Master Mixtures.

The recommend maximum freeze-thaw cycle is five cycles.

Additional Reagents and Equipment Required but Not Supplied

- 1) Compatible PCR instruments: AB 7500 / 7500 FAST, QuantStudio 5, CFX96, CFX384 Real-time system.
- 2) DNA Extraction kit:
(Recommended : QIAamp DSP DNA Blood Mini Kit for Whole blood, QIAamp MinElute ccfDNA kit for plasma sample)
- 3) Spectrophotometer for measuring DNA concentration.
- 4) Mini centrifuge with rotor for centrifuge tubes.
- 5) Mini centrifuge with rotor for PCR tubes.
- 6) Nuclease-free centrifuge tubes.
- 7) Nuclease-free PCR tubes and caps.
- 8) Adjustable pipettors and filtered pipette tips for handling DNA.
- 9) Tube racks.
- 10) Disposable powder-free gloves.
- 11) Sterile, nuclease-free water.

Precautions and Handling Requirements

For research use only.

Precautions

- DO NOT use for *in-vitro* Diagnostic (IVD).
- Please read the instruction carefully and become familiar with all components of the kit prior to use, and strictly follow the instruction during operation.
- Please check the compatible real-time PCR instruments prior to use.
- DO NOT use the kit or any kit component after their expiry date.
- DO NOT use any other reagents from different lots in the tests.
- DO NOT use any other reagent in the other test kits.

Safety Information

- Handle all specimens and components of the kit as potentially infectious material using safe laboratory procedures.
- Only trained professionals can use this kit. Please wear suitable lab coat and disposable gloves while handling the reagents.
- Avoid skin, eyes and mucous membranes contact with the chemicals. In case of contact, flush with water immediately.
- DO NOT pipet by mouth.

Decontamination and Disposal

- The kit contains positive control; strictly distinguish the positive control from other reagents to avoid contamination which may cause false positive.
- PCR amplification is extremely sensitive to cross-contamination. The flow of tubes, racks, pipets and other materials used should be from pre-amplification to post-amplification, and never backwards.
- Gloves should be worn and changed frequently when handling samples and reagents to prevent contamination.
- Using separate, dedicated pipettes and filtered pipette tips when handling samples and reagents to prevent

exogenous DNA contamination to the reagents.

- Please pack the post-amplification tubes with two disposable gloves and discard properly. DO NOT open the post-amplification PCR tubes.
- All disposable materials are for one time use. DO NOT reuse.
- The unused reagents, used kit, and waste must be disposed of properly.

Cleaning

- After the experiment, wipe down the work area, spray down the pipettes and equipment with 75% ethanol or 10% hypochlorous acid solution.

Instrument Setup

- Setup the reaction volume as 20 μ L.
- For AB 7500 / 7500 FAST, QuantStudio 5, CFX96, CFX384 Real-time system, please set up as follows: Reporter Dye: FAM, CY5; Quencher Dye: NFQ-MGB; Passive Reference dye: ROX.
- Refer to the real-time PCR instrument operator's manual for detailed instructions.
- We recommend that all PCR instruments in use should be conducted fluorescence calibration once a year.

Assay Procedure

1. DNA Extraction

The specimen material must be human genomic DNA extracted from whole blood samples or circulating DNA extracted from the plasma sample. DNA extraction reagents are not included in the kit. Before DNA extraction, it's essential to use standard pathology methodology to ensure tumor sample quality. Carry out the DNA extraction according to the instructions of DNA extraction kit.

The OD₂₆₀/OD₂₈₀ value of extracted DNA from whole blood sample or plasma sample should be between 1.7 ~ 2.1 (measured using the spectrophotometer). The storage conditions for extracted DNA from whole blood and plasma used for PCR amplification is shown in Table 3. And the circulating DNA isolated from plasma should be used directly without dilution.

Table 3. Recommended Conditions for Sample and Extracted DNA

	$\leq -70^{\circ}\text{C}$	$-15^{\circ}\text{C to } -25^{\circ}\text{C}$	$2^{\circ}\text{C to } 8^{\circ}\text{C}$	$15^{\circ}\text{C to } 30^{\circ}\text{C}$
whole blood sample			Up to 3 days	Up to 1 days
Plasma	Up to 12 months		Up to 3 days	
Extracted DNA		Up to 2 freeze thaws over 60 days	Up to 21 days	Up to 7 days

Note:

- The whole blood sample should be handled and stored properly, and the storage time should preferably be less than 3 days.
- The whole blood sample should be derived from EDTA anti-coagulated whole blood samples. The plasma sample should be derived from EDTA anti-coagulated peripheral whole blood. The recommended volume of whole blood is no less than 10 mL. The PAXgene Blood DNA Tube allows collection of 10mL whole blood into a closed

evacuated system.

- *The extracted DNA should be used immediately, if not, it should be stored at -20±5°C for no more than 2 months.*
- *Before detection, dilute the extracted DNA with 1×TE buffer (pH 8.0) to designated concentration. We recommend using at least 5 µL DNA for 10 times dilution, to ensure the validity of final concentration.*

2. Mutation Detection

- 1) Take the **JAK2 V617F Master Mixture**, **JAK2 V617F Positive Control**, **Nuclease-free water** and **ADPS smart DNA polymerase** out of the kit from the freezer.
- 2) Thaw the **JAK2 V617F Master Mixture**, **JAK2 V617F Positive Control** and **Nuclease-free water** at room temperature. When the reagents are completely thawed, vortex each tubes for 3 sec and centrifuge briefly to collect all liquid at the bottom of the tube.
- 3) Gently tap the **ADPS™ smart DNA polymerase** and centrifuge briefly prior to use.
- 4) Prepare sufficient **JAK2 V617F Reaction Mixture** containing **ADPS™ smart DNA polymerase**, **JAK2 V617F Master Mixture** in separate sterile centrifuge tube respectively according to the ratio in Table 4-1, 4-2. Mix Reaction Master Mix thoroughly by vortexing for 3 sec and centrifuge briefly.

Table 4-1. JAK2 V617F Reaction Master Mixture for Mutation Detection

Content	Volume per test (µL)
JAK2 V617F Master Mixture	10.0
ADPS™ smart DNA polymerase	0.5
Sample	2.0~9.5
Nuclease-free water	0~7.5
Total volume	20.0

Table 4-2. JAK2 V617F Reaction Master Mixture for Positive control

Content	Volume per test (µL)
JAK2 V617F Master Mixture	10.0
ADPS™ smart DNA polymerase	0.5
JAK2 V617F Positive Control	5.0
Nuclease-free water	4.5
Total volume	20.0

Note:

- *Prepare a positive control experiment at the very end to prevent cross-contamination.*
- *Every PCR run must contain at least one PC (Positive control) and one NTC (No template control).*
- *The prepared mixtures should be used immediately, avoid prolonged storage.*
- *Due to the viscosity of the enzyme mix, pipet slowly to ensure all mix is completely dispensed from the tip.*
- *Pipet enzyme mix by placing the pipet tip just under the liquid surface to avoid the tip being coated in excess enzyme.*

- 5) Prepare two PCR tubes for each sample: Dispense 10.0 µL of JAK2 V617F Reaction Mixture to each PCR tube respectively. Then add 5.0 µL of each sample DNA to each sample tube and cap the PCR tubes.
- 6) Prepare one PCR tube for PC: Dispense 10.0 µL of JAK2 V617F Reaction Mixture to one PCR tube. Then add PC (5 µL) to sample tube and cap the PCR tube.
- 7) Add nuclease-free water to all PCR tubes upto 20.0 µL.
- 8) Briefly centrifuge the PCR strips to collect all liquid at the bottom of each PCR tube.
- 9) Place the PCR strip tubes into the real-time PCR instrument. A recommended plate layout is shown in Table 5.

Table 5. PCR Plate Layout (Example)

Well	1	2	3	4	5	6	7	8	9	10	11	12
A	NTC	Positive Control										
B	Sample 1	Sample 8										
C	Sample 2	Sample 9										
D	Sample 3	Sample 10										
E	Sample 4	Sample 11										
F	Sample 5	Sample 12										
G	Sample 6	Sample 13										
H	Sample 7	Sample 14										

- 10) Setup the PCR Protocol using the cycling parameters in Table 6.

Table 6. Cycling Parameters

Stage	Cycles	Temperature	Time	Data collection
1	1	95°C	5 min	-
		95°C	10 sec	-
2	50	58°C	30 sec	FAM/Quasar670 (CY5)
		72°C	10 sec	-

- 11) Run the PCR run immediately.
- 12) When the PCR operation finished, analyze the data according to the “Results Interpretation” procedures

3. Results Interpretation

Before mutation data analysis, the following items should be checked:

- 1) Set the threshold to 10% of the maximum fluorescence value of the Positive control (JAK2 V617F).
- 2) For Positive control: FAM Ct value range of the samples should be 24 <Ct < 28.0. If not, the data is INVALID. The sample should be retested.
- 3) For the Internal control assay for each sample: Ct values of the samples should be 23 <Ct < 26 (Based on copy number of total JAK2 gene). If Ct value of IC < 23, this indicates that there are too many template in reaction. In this case, it would be better to dilute samples and should be retested. If Ct value of IC > 26, this indicates there are insufficient DNA or PCR inhibitors. The sample should be retested with increased or re-extracted DNA.

Analyze the mutation assay for each sample:

- 1) Record the FAM/Quasar670 Ct values of the samples.
- 2) Calculate the ΔCt value for each well:
 $\Delta Ct \text{ value} = \text{Sample Ct value (FAM)} - \text{Positive control Ct value (FAM)}$
- 3) ΔCt Cut-off value: $\Delta Ct < 19.0$
- 4) Result interpretation for each tube according to the Cut-off ΔCt value in Table 7.
 - a) If the ΔCt value is $<$ the Cut-off ΔCt value, the sample is determined as positive.
 - b) If there are no FAM signal and Ct values of sample is $44 < Ct$, the sample is determined as negative.
 - c) If the ΔCt value is \geq the Cut-off ΔCt value, the sample is determined as negative or under the LoD (limit of Detection) of the kit.

Table 7. Result Determination

JAK2 V617F Master Mixture			Determination
	FAM	Quasar670 (Cy5)	
Positive Control Ct value	$Ct < 24, 28.0 < Ct$	—	Invalid
Internal Control Ct value	—	$Ct < 23, 26.0 < Ct$	Invalid
Sample Ct value	$Ct < 23.5$		Invalid
	$23.5 \leq Ct < 45.0$		Positive
	$45 \leq Ct$		Negative
	No signal		Negative
ΔCt Cut-off value	$\Delta Ct < 19.0$		Positive

Performance Characteristics

The performance characteristics of this kit were validated on AB 7500 Fast.

1. Limit of Detection

The Low Limit of Detection (LLoD) of the kit for each mutation is shown in Table 8.

Table 8. LoD for each JAK2 V617F mutation

Codon	Mutation	Base Change	COSMIC ID	LoD (copy)	Detection sensitivity (%)
617	V600E1	G1849T	12600	3	0.01

2. Cross-reactivity

The cross reaction among the mutant sequences targeted by this kit, the cross reaction with other homologous mutant nucleotide sequence, the cross reaction with wild-type genomic DNA (DNA concentrations are 99 ng/reaction) were evaluated, the results shown no cross reactions.

3. Interference factor

1 interference substance: To evaluate the potential interference in interference substances on the performance of the ADPS™ JAK2 V617F Mutation Test kit (RUO). The levels of potential interfering substance were equal to the levels recommended to be tested by the Clinical and Laboratory Standards Institute (CLSI) EP7-A2. It is confirmed that the potential maximum concentrations: 30 copies/5 µ would not interfere with the test result.

4. Precision

3 precision controls: negative control, weak positive control (with 0.01% mutant content) and strong positive control (with 1% mutant content) were used in the validation. 3 batches of the kits were tested with the precision controls by 2 operators twice a day for 20 days on different PCR instruments. The Ct values were calculated, the CV values were all within 5%.

Limitations

- 1) The kit is to be used only by personnel specially trained with PCR techniques.
- 2) The 0.01 % of detection sensitivity is guaranteed only when input DNA is used with maximum of 30,000 copies.
- 3) The results should never be used for *in vitro* diagnosis (IVD).
- 4) The kit has been validated for use with human genomic DNA extracted from whole blood sample and circulating DNA extracted from plasma sample.
- 5) Reliable results are dependent on proper specimen collection, processing, transport, and storage.
- 6) The sample containing degraded DNA may affect the ability of the test to detect *JAK2 V617F* mutation.
- 7) This kit can only assess the *JAK2 V617F* mutation status and detect *JAK2 V617F* mutation indicated above.
- 8) Samples with negative result (No mutation detected) may harbor another *JAK2* mutations not detected by this assay.









References

- 1) Vardiman JW, Thiele J, Arber DA, et al. The 2008 revision of the World Health Organization (WHO) classification of myeloid neoplasms and acute leukemia: rationale and important changes. *Blood* 114 (5): 937-51, 2009
- 2) Arber DA, Orazi A, Hasserjian R, et al.: The 2016 revision to the World Health Organization classification of myeloid neoplasms and acute leukemia. *Blood* 127 (20): 2391-405, 2016
- 3) Savona MR, Malcovati L, Komrokji R, et al. An international consortium proposal of uniform response criteria for myelodysplastic/myeloproliferative neoplasms (MDS/MPN) in adults. *Blood* 125 (12): 1857-65, 2015.
- 4) James, C., Ugo, V., Le Couedic, J.P., et al. A unique clonal *JAK2* mutation leading to constitutive signaling causes polycythaemia vera. *Nature* 2005. 434, 1144.
- 5) Tan AY et al. A Simple, Rapid, and Sensitive Method for the Detection of the *JAK2 V617F* Mutation. *Am J Clin Pathol.* 2007 Jun;127(6):977-81.
- 6) Rane SG, Reddy EP. Janus kinases: components of multiple signaling pathways. *Oncogene.* 2000;19:5662-5679.
- 7) Tefferi A, Gilliland D. The *JAK2V617F* tyrosine kinase mutation in myeloproliferative disorders: status report and immediate implications for disease classification and diagnosis. *Mayo Clin Proc.* 2005;80:947-958.
- 8) Klampfl T, Gisslinger H, et al. Somatic mutations of calreticulin in myeloproliferative neoplasms. *N Engl J Med.* 2013 Dec 19;369(25):2379-2390
- 9) Tan, A. et. al. (2007 July 23). A Simple, Rapid, and Sensitive Method for the Detection of the *JAK2 V617F* Mutation. Medscape from *American Journal of Clinical Pathology* [On-line information]. Available online F701-19(0)

at <http://www.medscape.com/viewarticle/558906> through <http://www.medscape.com>. Accessed July 2009.

- 10) Tefferi et al. Polycythemia vera and essential thrombocythemia: 2017 update on diagnosis, risk-stratification, and management. *Am J Hematol.* 2017 Jan;92(1):94-108.
- 11) Schischlik et al. Mutations in myeloproliferative neoplasms- their significance and clinical use. *Expert Rev Hematol.* 2017 Nov;10(11):961-973.
- 12) Rumi et al. Diagnosis, risk stratification, and response evaluation in classical myeloproliferative neoplasms. *Blood.* 2017;129(6):680-692.
- 13) Wang et al. A portable microfluidic platform for rapid molecular diagnostic testing of patients with myeloproliferative neoplasms. *Nature. Scientific Reports.* 7: 8596. 2017.

Symbol Information

Symbol	Title of symbol	Symbol	Title of symbol
	Temperature limit		Caution
	Use-by date		Consult instructions for use
	Batch code		
	Catalogue number		
	Manufacturer (Name&Address)		
	Contains sufficient for <n> tests		